

## Filter Connectors for 25mm & 50mm diameter filters.

### 25mm Diameter Filter Connectors:

In-Let	Out-Let
Female Luer Lock	Male Slip Luer
Female Luer Lock	Male Luer Lock
Female Luer Lock	Female Luer Lock
Female Luer Lock	3/16" ~ 1/4" Stepped Barb
3/16" ~ 1/4" Stepped Barb	3/16" ~ 1/4" Stepped Barb
3/16" ~ 1/4" Stepped Barb	Male Slip Luer
3/16" ~ 1/4" Stepped Barb	Male Luer Lock
3/16" ~ 1/4" Stepped Barb	Female Luer Lock
Male Luer Lock	Male Luer Lock
Male Luer Lock	Female Luer Lock
Male Luer Lock	Male Slip Luer
Male Luer Lock	3/16" ~ 1/4" Stepped Barb
Male Slip Luer	Male Slip Luer
Male Slip Luer	3/16" ~ 1/4" Stepped Barb
Male Slip Luer	Male Luer Lock
Male Slip Luer	Female Luer Lock

### 50mm Diameter Filter Connectors:

In-Let	Out-Let
Female Luer Lock	3/16" ~ 1/4" Stepped Barb
Female Luer Lock	Female Luer Lock
Female Luer Lock	Male Luer Lock ~ Finless
Female Luer Lock	1/4" - 3/8" Stepped Barb
Female Luer Lock	Flush Connector
Female Luer Lock	1/8" British Standard Pipe Thread
Female Luer Lock	1/8" National Pipe Thread
Male Luer Lock ~ Finless	Male Luer Lock ~ Finless
Male Luer Lock ~ Finless	1/4" - 3/8" Stepped Barb
Male Luer Lock ~ Finless	3/16" ~ 1/4" Stepped Barb
Male Luer Lock ~ Finless	Flush Connector
Male Luer Lock ~ Finless	1/8" British Standard Pipe Thread
Male Luer Lock ~ Finless	3/8" ~ 1/2" Stepped Barb ~ Finless
1/4" - 3/8" Stepped Barb	1/4" - 3/8" Stepped Barb
1/4" - 3/8" Stepped Barb	3/16" ~ 1/4" Stepped Barb
1/4" - 3/8" Stepped Barb	Female Luer Lock
1/4" - 3/8" Stepped Barb	Male Luer Lock ~ Finless
3/16" ~ 1/4" Stepped Barb	3/16" ~ 1/4" Stepped Barb
1/8" National Pipe Thread	Flush Connector
1/8" National Pipe Thread + Jaco Connector	Flush Connector
1/8" British Standard Pipe Thread	1/8" British Standard Pipe Thread
1/8" British Standard Pipe Thread	Flush
5/16 Unified National Fine Thread + Ferrell Nut	5/16 Unified National Fine Thread + Ferrell Nut
Flush	3/8" ~ 1/2" Stepped Barb ~ Finless
1/4" - 3/8" Stepped Barb	3/8" ~ 1/2" Stepped Barb ~ Finless

All connectors are 'Finned' unless otherwise stated. Fins are used to add strength to the filter & connector.

We are also able to produce 'Coloured' filters upon request. Please note there is a Minimum Order Quantity on Coloured Filters.

If you cannot see your desired combination of connectors above, please contact us for details.

*(Helapet reserve the right to modify this list at any time, and by providing advice on the workable combinations of connectors)*

**HELAPET**

Helapet Limited | Circle Business Centre | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5DD | United Kingdom  
Telephone +44 (0) 1582 501980 | Fax +44 (0) 1582 501981 | e-mail sales@helapet.co.uk | www.helapet.co.uk