

HELAPET

VARIDISK™


- ❖ Upstream Support
Choice of Inlet Connector
to suit you application
Fins for extra pressure
resistance
Optional upstream Vent

- ❖ Filter Media
Your choice of media

- ❖ Helapet Usage Ring
Unique optional
time/date/usage ring

- ❖ Sealing Ring
Rugged construction

- ❖ Downstream Support
Filter media sealed to
eliminate bypass
Choice of Outlet Connector



Custom Filter Devices

A custom 50mm device built
to your specification.

FILTRATION

"Filtration by innovation"